


PHYSION[®]



N-DO[®]

THE NEW ENDOINJECTOR SYSTEM
for day-surgery endoscopic procedures



N-DO[®] Injector

- Curving needle with slider and safety indicator
- The distal Nitinol needle is adjustable in length up to 16 mm.



- Due to its geometrical characteristics, **N-DO[®]** needle is always visible which guarantess **safety, effectiveness and control**. **N-DO[®]** needle allows **up to 30 injections** in one patient.


- Thanks to its easy handling, **N-DO[®]** needle **allows to perform a precise operation** in different areas of the low urinary system, otherwise difficult to treat (e.g. bladder wall, bladder neck).

N-DO[®] is a catheter curved needle with **slider and safety indicator**

to be used with a cystoscope that is **adjustable by length** in increments of 2mm,

allowing selected drug solutions to be injected at the **desired depth**

in a **safe and controlled manner**.



INDICATIONS

- Local anaesthesia:
Minor TURBT procedures (transurethral resection of bladder tumour)
TURP - (transurethral resection of the prostate)
- Bladder biopsies
- Injections of drugs and chemotherapies directly into the bladder
- Toxic Botulinum injections



N-DO® Injector for use with rigid cystoscope
45 cm



N-DO® Injector for use with flexible cystoscope
70 cm



N-DO[®] Injector
provides a controlled,
easy and quick
method to perform
day surgery
endoscopic
procedures.

PHYSION[®]

Physion Srl

Via Fogazzaro 4/B 41037 Mirandola (MO) - Italy

Phone: +39 0535 611082 Fax: +39 0535 610420

email: info@physion.it

physion.it