


UVENTATM Urethral Stent

for benign & malignant urethral stricture



**SHORT-TERM USE,
LONG-LASTING EFFECT**

HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Gyeonggi-do, Korea 10355

E: lea@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

UVENTA™ Urethral Stent

Indication



- Benign and malignant urethral obstruction
- at Bulbous and Penile urethral sites

Features

- Fully-covered with Silicone
- Self-expandable
- Anti-migration barbs
- Radiopaque markers
- Retrieval Suture

Delivery system


- Pre-mounted on the delivery system



- **Anti-migration barbs**
to prevent stent migration



- **Retrieval suture**
for easy removal of stent



Ordering Information

Specification of UVENTA™ Urethral Stent

Code	Stent		Delivery		Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
US1004F	10	4	10	50	US1404F	14	4	10	50
US1005F		5			US1405F		5		
US1006F		6			US1406F		6		
US1008F		8			US1408F		8		
US1010F		10			US1410F		10		
US1204F	12	4			US1604F	16	4		
US1205F		5			US1605F		5		
US1206F		6			US1606F		6		
US1208F		8			US1608F		8		
US1210F		10			US1610F		10		