




Soft lift

The Stress Urinary Incontinence solution for all your patients




1 kit : **4 surgical techniques**

Precise sling adjustment

Optimised patient comfort

SOFT LIFT : THE STRESS URINARY INCONTINENCE SOLUTION FOR ALL YOUR PATIENTS

> KEY POINTS



Suprapubic approach &
Retropubic approach

1 kit : **4 surgical techniques**

- TVT and TOT approaches



TOT In-Out approach &
TOT Out-In approach

Optimised patient comfort

- Light weight polypropylene
- Smooth edges

> REFERENCES AND SIZES

Reference	Description and size (cm)
GCB INSOFT U	Single use kit SOFT LIFT

> EACH SINGLE-USE KIT INCLUDES:

1 Implant	2 TOT Ancillaries	1 Curved Ancillary	1 Guide
 (1 cm x 45 cm)			

> TECHNICAL CHARACTERISTICS

- Knitted polypropylene monofilament
- Macroporous
- Light weight ($64 \pm 7 \text{ g/m}^2$)
- Low elasticity
- Interstice ($0,003 \text{ mm}^2$)
- Pore size ($0,8 \times 0,9 \text{ mm}$)

Intramesh SOFT LIFT is a class IIb medical device manufactured by Cousin Biotech s.a.s. The management system of COUSIN BIOTECH s.a.s has been certified as meeting the requirements of ISO 13485. Please read instructions for use carefully. Reference: FPGSLLGB02 - 19/01/17.

Non contractual pictures and texts. Specifications likely to be modified without notice.
Cousin Biotech S.A.S capital : 340 656 € - 398 460 261 RCS Lille - N°TVA FR 34 398 460 261

CE
0120

