

Testi10™ testicular prosthesis

Testi10™ is designed with high quality silicone. The prosthesis is shaped and designed with a soft texture to offer patients missing one or both testicles a more natural like feeling.



Saline Filled
Testicular Prosthesis

- ♂ More natural feeling designs
- ♂ Aims for higher patient satisfaction
- ♂ Available in five sizes to better meet the patient's anatomical needs



TestiSF-XL



TestiSF-L



TestiSF-M



TestiSF-S



TestiSF-XS



Testi10™ Saline-filled Testicular Prosthesis		
Prosthesis Product	Dimension A	Dimension B
TESTISF-XS	2.7 cm	2.1 cm
TESTISF-S	3.2 cm	2.6 cm
TESTISF-M	4.1 cm	2.9 cm
TESTISF-L	4.7 cm	3.2 cm
TESTISF-XL	5.1 cm	3.2 cm

Testicular Prosthesis

Testi10™

Designed using high quality silicone, to help patients experience a complete and natural feeling.

Testi10™ is specifically used for the replacement of one, or two testicles due to monorchism, orchiectomy, congenital anorchism or similar procedures.

Saline-filled model is supplied as empty and gives the surgeon the opportunity to create a more natural feel to the implant

Testi10S™ At A Glance

Dedicated Injection Port

Dedicated injection ports for easy filling of the device during the surgery.

Diverse Sizes

Available in 5 sizes between XS to XL to offer the best fit to patient's anatomy.



Easy to Fill

Dedicated injection ports for easy filling of the device during the surgery



Diverse Sizes

Available in 5 sizes between XS to XL to offer the best fit to patient's anatomy.



Mimics the Natural

Mimics the Natural The prosthesis is designed to replicate the structure and feel of the natural testicle.



Robina
Group Companies

www.urologic.ir

(۰۲۱) ۴۴۰ ۲۵ ۴۹۵, ۴۴۰ ۲۵ ۲۰۷

بزرگراه ستاری شمال، خیابان پیامبر مرکزی، خیابان مطهری، پلاک ۵، واحد ۷

Rigicon

American Urological Solutions

NY 11779 United States of America