

SEESHEEN

Cy


 Image



**CY** series  
Flexible Video Cystoscope

## Suction function

One-button control to maintain clear view.  
Unique anti-drop design to prevent the suction valve from falling off.



\*Patent NO.: ZL202022141662.4

## Angle lock design

Lock or release the angle of the distal end for accurate diagnosis and treatment.



\*Patent NO.: ZL202022141661.X

## Multi-function buttons

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release



## Tapered tip


Tapered tip design for easier insertion.




### Product Parameter

Model	Insertion Tube	Working Channel	Working Length	Bending (Up/Down)	Rotation (Left/Right)	Field of View	Depth of field	Suction function
CY-1356	15.6Fr	6.6Fr	380mm	220°/130°	120°	120°	3-100mm	Yes
CY-1355	15.6Fr	6.6Fr	380mm	220°/130°	120°	120°	3-100mm	No


### Image Processor




◀ HV-3101  
Image Processor




10.1"  
Touch screen




1280\*800  
HD Resolution




HDMI+ AV  
output



Detachable & rechargeable battery  
supporting 4 working-hour



Detachable  
SD card



Portable  
design