

SEESHEEN

UR



UR series

Flexible Video
Ureterorenoscope


Secondary bending

Bending U/D 275°
Secondary bending to improve operation flexibility



Rotation

Improve operation flexibility and accuracy



*Patent NO.: ZL201621392021

Tapered tip

Tapered tip design for easier insertion



Multi-function buttons

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release



Product Parameter

Model	Insertion Tube	Working Channel	Working Length	Bending (Up/Down)	Rotation (Left/Right)	Field of View	Depth of Field
UR-1328	8.5Fr	3.6Fr	670mm	275°/275°	120°	120°	3-50mm

Image Processor



10.1" Touch screen



1280*800 HD Resolution



HDMI+ AV output



Rechargeable battery supporting 4 working-hour



Detachable 64G SD card



Portable design

