


Disposable Flexible Uretroscope

Twin light guides deliver a clear and bright image.



OD.: 3.0mm (9Fr.)

OD.: 2.5mm (7.5Fr.)

Reduces the risk of high kidney pressure

- Slim 7.5Fr
- 275° deflection
- 400x400 pixels
- Oval Shape




7.5 Fr.

3.6 Fr. (1.2mm) Working Channel

Single use
Super image quality
Variety of models





Irrigation

Instrument

Cystoscope

less weight and ergonomic design
make the scope procedure more
easier to control.



Insertion Dia.

UV-100-E 14Fr.(4.8mm)
UV-100T-E 16 Fr. (5.3mm)



Working channel

UV-100-E 6.6Fr.(2.2mm)
UV-100T-E U7.8 Fr.(2.6mm)

- Power**

Voltage 100-240 VAC
Frequency 50/60 HZ
Rared input 45 VA

- Illumination**

Lamp 2W LED
Color temperature Lamp life 30000 hrs
3000k- 7000k

- Function**

White balance (+, -), Illumination control, Enhance
(0,1,2), Freeze, Zoom (x1,x2,x3),
Snap, Record (USB Storage)

- Output**

DVI*1, S-Video*1, CVBS*1

- Storage**

USB (FAT32)

- Dimension**

390mm(W)* 270mm(H)* 200mm(D)

- Weight**

7.2 kg

- Screen**

12 inch

Imaging processor PL-2100



160k high resolution image

The video processor
incorporates both freeze and
video capture to USB.

SEEGEN

Shanghai SeeGen Photoelectric Technology Co.,Ltd.



Robina
Group Companies

Telefax: (021) 440 25 207

Unit 07, No 05, Motahari St, Central Payambar St, Tehran, Iran

www.urologic.ir

